Omnidouble 15

**Additionsvernetzendes 1 : 1 Dubliersilikon
(DIN EN ISO 14356, Typ 2, irreversible Dubliermassen)**

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Omnidouble 15

Addition curing 1:1 duplicating silicone

(DIN EN ISO 14356, Typ 2, irreversible duplicating material)

Technical data:	Mixing ratio A (white) : B (blue):	1 : 1
	Mixing under vacuum:	40 seconds
	Automatic measuring:	in a dispenser with a mixing nozzle*
	Working time at 23°C:	> 5 minutes
	Curing time at 23°C:	30 - 45 minutes
	Tensile strength:	> 1.8 MPa
	Elongation at break::	> 220 %
	Tear strength:	approx. 3.0 N/mm
	Shore A hardness after 60 minutes	approx. 15

Application: Omnidouble 15 is suitable for all duplicating flasking systems. Models to be duplicated should be at room temperature and should not be soaked. Any dirt on the surface is removed by spraying with an alcohol based neutralizer and by blow-drying carefully. Materials, such as gypsum hardener, grease, oils and some separators, can produce a negative reaction with Omnidouble 15 (refer to “Incompatible materials” below). When selecting the flask size, ensure there is a gap of at least 5mm between the model and the flask so that the Omnidouble 15 can flow round the model. Mix Omnidouble 15 under vacuum and pour in a thin flow until the model is completely covered. The stabilising or base plate is then placed in position. Omnidouble 15 can be measured using a 1:1 dispenser because both components have almost the same viscosity.

* If a 1:1 dispenser is used for measuring, curing should be carried out under pressure at 2.5 – 3.0 bar. Use the duplicating cross instead of the fixation device to secure the stabilising plate when curing under pressure.

Chemical characterisation: Polymerised dimethyl polysiloxane filler mixture

Dangerous reactions: No known dangerous reactions, if stored and used properly.

Toxicity: None known

Directions: Store containers tightly sealed in a well-ventilated place. Keep containers dry. Keep away from water, alkalis, strong acids and oxidising agents.

Ecological information: Present experience suggests that there is no detrimental effect on the environment, if it is used properly and current regulations are observed. It is non-biodegradable and insoluble in water when cured. It can be easily separated from water by filtering.

Please note: Cured silicone is not harmful to the environment and can be disposed of along with domestic waste, as long as this does not contravene local bye-laws. Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. OMNIDENT products are subjected to strict quality controls. We reserve the right to make technical changes.

Incompatible materials:

The following materials can produce a negative reaction with addition-curing silicones: soft polyvinyl chloride, polysulphide MIL-S-8516, Mystik[®] 6207 adhesive tape, Mystik[®] 6215 adhesive tape, Scotch[®] cellophane strips, Scotch[®] 360 adhesive tape, Permacel[®] masking tape, Pliobond[®] adhesive, Neoprene rubber, Buna N rubber, Vitron[®]A rubber, colophony-based soldering flux, acid-based soldering flux, disinfectant, ISOLIT, cyanoacrylate glue, amino-hardened epoxy, Humiseal[®] 6215 adhesive tape, vinyl insulating tape, latex vacuum hoses, GRS rubber, natural gum, certain types of silicone and RTV rubber, SUPER-SEP, milling oils, sulphur compounds such as thiols, sulphides, sulphates, sulphites, thiourea and nitrogen compounds such as amines, amides, imides, acids.

If there is any doubt, we recommend testing the material before use.

Rev. 2204_02 / 2022-04